

## Extreme mobility: Low-temperature NMR probes highly diffusive Li<sup>+</sup> ions in argyrodite-type Li<sub>6</sub>PSe<sub>5</sub>Cl and Li<sub>6</sub>PS<sub>5</sub>Br

Viktor Epp<sup>1\*</sup>, Özgül Gün<sup>2</sup>, Hans-Jörg Deiseroth<sup>2</sup>, Martin Wilkening<sup>1</sup>

<sup>1</sup>Christian Doppler Laboratory for Lithium Batteries, Institute for Chemistry and Technology of Materials, Graz University of Technology, Graz, Austria

<sup>2</sup>Anorganische Chemie, University of Siegen, Siegen, Germany

\*viktor.epp@tugraz.at

Lithium argyrodites represent a new class of crystalline solids exhibiting a very high Li<sup>+</sup> conductivity and a negligible electronic conductivity [1, 2]. This makes them ideally suited to act as electrolytes in highly stable and reliable all-solid-state lithium-ion batteries. In general, the use of a solid electrolyte inherently resolves a large majority of the safety issues associated with this technology.

The accurate characterization of Li dynamics in crystalline solids helps identify the most promising candidates so that subsequent research can focus on the improvement of their overall properties for battery applications. Here, Li diffusion parameters in a series of Li-argyrodites have been comprehensively studied by using a range of complementary time-domain NMR techniques. In particular, these include the acquisition of diffusion-induced NMR spin-lattice relaxation (SLR) rates in both the laboratory ( $R_1$ ) and in the rotating ( $R_{1\rho}$ ) frame of reference [3, 4].




Figure 1: Arrhenius plot of the  ${}^7\text{Li}$  SLR NMR rates recorded with the spin-lock technique in the rotating frame of reference at frequencies of 14 kHz and 40 kHz, respectively; taken from Ref. [4]. The parameter  $\beta$  ( $1 < \beta \leq 2$ ) indicates the deviation from ideal BPP-type behaviour [5] which would be reflected by  $\beta = 2$ .

Most importantly, for the Cl- and Br-containing compounds,  $\text{Li}_6\text{PSe}_5\text{Cl}$  and  $\text{Li}_6\text{PS}_5\text{Br}$ , extremely fast Li exchange processes are found. This manifests in diffusion-induced (laboratory frame) SLR NMR rate peaks showing up at temperatures as low as 260 K while the corresponding rate peaks measured in the rotating frame of reference are shifted towards much lower temperatures (see figure and Ref. [4]). From a quantitative point of view, at ambient temperature the Li jump rate in polycrystalline  $\text{Li}_6\text{PS}_5\text{Br}$  is of the order of  $10^9 \text{ s}^{-1}$  which corresponds to an Li ion conductivity in the order of  $10^{-2} \text{ S/cm}$ , thus, indicating ultrafast translational motion of the charge carriers.

## References

- [1] H.-J. Deiseroth, S.-T. Kong, H. Eckert, J. Vannahme, C. Reiner, T. Zaiß, M. Schlosser: *Li<sub>6</sub>PS<sub>5</sub>X: A class of crystalline Li-rich solids with an unusually high Li<sup>+</sup> mobility*. *Angew. Chem. Int. Ed.* **47**, 755–758 (2008)
- [2] S.-T. Kong, Ö. Gün, B. Koch, H.-J. Deiseroth, H. Eckert, C. Reiner: *Structural characterisation of the Li argyrodites Li<sub>7</sub>PS<sub>6</sub> and Li<sub>7</sub>PSe<sub>6</sub> and their solid solutions: Quantification of site preferences by MAS-NMR spectroscopy*. *Chem. Eur. J.* **16**, 5138–5147 (2010)
- [3] V. Epp, Ö. Gün, H.-J. Deiseroth, M. Wilkening: *Long-range Li<sup>+</sup> dynamics in the lithium argyrodite Li<sub>7</sub>PSe<sub>6</sub> as probed by rotating-frame spin–lattice relaxation NMR*. *Phys. Chem. Chem. Phys.* **15**, 7123–7132 (2013)
- [4] V. Epp, Ö. Gün, H.-J. Deiseroth, M. Wilkening: *Highly mobile ions: Low temperature NMR directly probes extremely fast Li<sup>+</sup> hopping in argyrodite-type Li<sub>6</sub>PSe<sub>5</sub>Br*. *J. Phys. Chem. Lett.* **4**, 2118–2123 (2013)
- [5] N. Bloembergen, E.M. Purcell, R.V. Pound: *Relaxation effects in nuclear magnetic resonance absorption*. *Phys. Rev.* **73**, 679–712 (1948)